

QA4

Surgery Trolley System



***The complete multi-purpose
mobile operating system***

QA4 *POWERED* functions

The QA4 Surgery System, both patient trolley and operating table, a total solution for modern day treatment centres.



- Electrically Powered Functions •
- Superior Surgical Access •
- Traversing Patient Platform •
- 100% 'C' Arm Access •
- Lateral tilt •
- 250kg Weight Limit •

MANUAL *functions* QA4

Choose either the powered or manual version to suit your needs.



- Mechanically Powered Functions
- Increased Surgical Access
- Supremely Lightweight & Manoeuvrable
 - 75% 'C' Arm Access
 - Low Cost Solution
 - 160kg Weight Limit



QA4 POWERED functions

With four remotely-activated powered functions, including a traversing top, the QA4 Powered trolley offers full 'C' arm coverage and superior surgical access together with a 250kg patient weight limit and exceptional stability.




TRANSPORT & RECOVERY

- Cotsides fitted for patient comfort and safety
- Unique retractable steering mechanism engaged for maximum manoeuvrability



• ± 12" Lateral tilt



• Unique traversing top allows for 100% 'C' arm access



GYNAECOLOGY

- Patient platform traversed to foot end for access
- Leg section removed for access
- Head section removed for anaesthetist access
- Fitted with optional powered leg supports



OPHTHALMICS

- Patient platform traversed head end for access
- Ophthalmic head section fitted



SHOULDER ARTHROSCOPY

- Patient platform traversed foot end for access
- Optional arthroscopy backrest fitted



Warranty subject to on-site servicing

COMMON FEATURES



Removable head section



Push button release and self locking cot sides



Push button rotating cot sides



Push button release and self locking leg section



Removable ultra lightweight leg section



K8 Pressure Care mattress



QA4 ACCESSORY trolley

The complete storage solution for general operating theatre and QA4 accessories



Warranty subject to on-site servicing

GYNAECOLOGY

- Lightweight leg section removed for access
- Head section removed for anaesthetist access
- Fitted with optional powered leg supports



ENT / DENTAL

- Patient positioned on standard articulating head section



- Low trolley height allows for easy patient mounting and dismounting



• 75% 'C' arm access

QA4 MANUAL functions



Excellent manoeuvrability and easy-to-use hand operated functions make the QA4 Manual trolley a cost effective option; a patient transport trolley with the features of an operating table.



TRANSPORT & RECOVERY

- Cotsides fitted for patient comfort and safety
- Unique retractable steering mechanism engaged for maximum manoeuvrability



SeeDOS Ltd

QA4 Surgery Trolley System



TROLLEY SPECIFICATIONS	QA4 POWERED Cat No. 21300	QA4 MANUAL Cat No. 21310
Trolley Features		
Powered Functions	✓	×
Lateral Tilt	✓	×
Head End Traverse	✓	×
Foot End Traverse	✓	×
K8 Pressure Care Mattress	✓	✓
Removable Ultra Lightweight Leg Section	✓	✓
Removable and Rotating Cotsides	✓	✓
Push Button Release and Self Locking Controls	✓	✓
5th Wheel Steering Mechanism	✓	✓
Advanced 4 Castor Braking System	✓	✓
'E' Size Oxygen Cylinder Trough (Side Loading)	✓	✓
Transfusion Pole with Mounting Position	✓	✓
150mm Castors	✓	✓
Height (MIN)	605mm	575mm
Height (MAX)	1005mm	935mm
Trolley Length	2100mm	2100mm
Side Rail Width	655mm	655mm
Cotside Width	770mm	770mm
Brake Width (Brakes Off)	835mm	835mm
Mattress Length	2040mm	2040mm
Mattress Width	600mm	600mm
Mattress Depth	75mm	75mm
Base Length	1190mm	1055mm
Trendelenberg	20°	20°
Reverse Trendelenberg	12°	-
Lateral Tilt	+12°/-12°	-
Backrest Articulation	0 - 65°	0 - 80°
Head Section Articulation	+30°/-30°	+30°/-30°
Head End Traverse	230mm	-
Foot End Traverse	230mm	-
'C' Arm Access	100%	75%
Max. Patient Weight Limits	250kg	160kg
Trolley Weight	160kg	125kg
Lightweight Leg Section Weight	6kg	6kg

Optional Extras

	Cat No.
QA4p Spare Transfusion Pole	21312
QA4m Spare Transfusion Pole	21313
Sculpted Mattress for Head Support	21315
Ophthalmic Head Support	21320
Removable Shoulder Arthroscopy Backrest	21340
Padded Cotside Covers - Grey	21350
Padded Cotside Covers - Child Print	21351
Articulating Removable Leg Section	21355
Removable Leg Section c/w Full Length Side Rail	21357
K8 Pressure Mattress Complete	21360
QA4 Accessory Trolley	21370
QA4 Accessory Trolley Spare Trough	21371
Spare Powered Trolley Hand Control	21380
Spare Powered Trolley Charging Cable	21385
CD Size Cylinder Support Bracket	53556



SeeDOS Ltd