

Intracavitary HDR Brachytherapy

CT / MRI Compatible 3 Channel Endometrial Applicator

Innovative Technology and Enhanced Clinical Advantages.
Facilitates CT and MRI based Brachytherapy Planning.

3 Channel Endometrial Applicator



Catalog # 0910

Shown: Varian - GammaMed® / VariSource® Applicator

Compatible With All of Today's Remote Afterloading Systems

The **Endometrial Applicator** is entirely constructed of medical grade Titanium.

Each tandem can be inserted separately and subsequently locked securely in place by means of only one lock-knob.

Lateral tandems can be locked in place with or without the central tandem.

The **Endometrial Applicator** is lightweight yet stable and sturdy.
