



www.seedos.com

Radiation Oncology Equipment and Services

Seed-Lock³

The Only Synthetic plugged needle for Brachytherapy

SeedLock³ is the better choice for pre-plugged needles. Its unique design means smoother loading, easier delivery, accurate placement, and reliable expulsion.

Protect your patient, protect yourself.

- Third generation pre-plugged needles utilize an advanced synthetic bio-absorbable micro-plug which provides superior retention and smoother expulsion.
- Consistent plug location allows for precise seed placement, taking full advantage of today's advanced planning and imaging modalities.
- New Accu-Trac™ needle tip design minimizes trauma while providing unequalled penetration and tracking.
- Ultra-enhanced echo tip improves imaging.
- Proprietary EZ Load™ technology virtually eliminates problems associated with seed/spacer loading and placement.
- Exclusive ESD clear hub material actively removes static from the spacers/seeds. Available in round or square hub.
- CM marks on cannula and .5cm marks on stylet configured for either standard or Seattle technique (reversed markings).
- Vented stylet provides an air passage to needle tip equalizing pressure to eliminate unwanted seed displacement from vacuum phenomena.

Seed-Lock³ is the only approved needle for use with Advanced Care Medical Vari-Strand™, Vari-Load™ and Real-Time Strands™

You have found the new standard of care.

www.wmedtech.com • Tel: 877-783-5463 • 203-262-4190
115 Hurley Rd. Building 3B • Oxford, CT, 06478 • Fax Orders: 800-789-0897

CE

0459 ISO 13485

Pre-plugged needles and pre-loaded needles are protected under US Patent No. 6,554,760.



