

Optional Shipping Case for QUASAR Beam Geometry Phantom or QUASAR Body Phantom

Designed to protect the QUASAR phantoms during repeated shipping and rough handling, the optional shipping case provides more protection than the storage and handling container that is provided with the phantoms.

Manufactured by Pelican Products of Torrance California (www.pelican.com) the shipping case is watertight, airtight, dustproof, and crushproof. Features include multiple fold down carrying handles, atmospheric pressure purge valve, and integral wheels. The shipping case is available from Modus, and approved distributors, with a custom designed closed cell foam insert to protect QUASAR anywhere you take it.



Case Specifications

Pelican Model 1620

Outside dimension: 24 13/16" L x 19 3/8" W x 13 7/8" D

Colour: Black



QUASAR Body Phantom
Model QBDP-case

QUASAR Beam Geometry Phantom
Model QBGP-case



Distributed by:

SeeDOS Ltd